

# Characterization and Potential Applications of Human iPSC-Derived Neural Stem Cells


Dutta-Passecker P<sup>1</sup>.; Gibbons G<sup>1</sup>.; Wessel-Carpenter N<sup>1</sup>.

<sup>1</sup>Axol Bioscience Ltd., UK. p.dutta@axolbio.com



We conducted a series of experimental procedures to examine the characteristics and potential application of iPSC-derived neural stem cells for neurobiological research.

## Introduction




- Adult cells can be reprogrammed into iPSCs using defined factors - Oct3/4, KLF4, Sox2 and c-Myc<sup>1-3</sup>.
- iPSCs can be differentiated into many cell types including human neural stem cells (hNSCs) and cerebral cortical neurons (hCCNs) from healthy donors and patients suffering from disease<sup>4-6</sup>.
- We characterized iPSC-derived hNSCs and their progeny to examine their suitability for neurobiology research.
- We have demonstrated that these cells can be used as a model for the study of human neuronal development.

## Results


### Transcriptome Analysis

Integration-free iPSCs were generated using an episomal vector and subsequently differentiated into hNSCs using Axol's proprietary method.



Experiments were performed using the Affymetrix GeneChip® Human Transcriptome Array 2.0 platform. Results were analysed using Affymetrix® Expression Console™ and Affymetrix® Transcriptome Analysis Console (TAC) 2.0 Software.

### Neural Cell Morphology




**Human Neural Stem Cells**  
Axol human neural stem cells (hNSCs) form neural rosettes and express markers typically observed in neural precursor cells as seen by immunocytochemistry.

**Human Cerebral Cortical Neurons**  
Differentiation of hNSCs generates human cerebral cortical neurons (hCCNs) that express neuronal markers observed using immunocytochemistry. These neurons increase in maturity over time in culture.

### 3D Culture


**Foxg1 - Nestin**  
On culturing hNSCs on top of the RAFT™ collagen matrix, (TAP Biosystems/Lonza) cell migration into the matrix was observed. The cells form the 3D structure of commonly seen neural rosettes. Outside of the rosette, there is a matrix of cells ordered in a non-uniform manner.

**Tbr1 - Tuj1**  
On culturing hCCNs on top of the collagen gel they formed a uniform static layer of cell bodies. Neurites projected out of these cells and grew downwards, creating a network of interconnected neurites.




### Neurite Outgrowth

Axol hNSCs yielded the highest neurite length and branch point values in comparison to rat primary cortical neurons (CTX), competitor iPSC-derived neurons, and N2a cells.



Axol hNSC neurite outgrowth was assessed by S. Lopez Alacantara & T. Dale, Essen Bioscience Ltd using the IncuCyte NeuroTrack platform.


### Whole Cell Patch Clamp




Whole cell patch clamp recordings were carried out by Ana González Rueda, Ole Paulsen Lab, University of Cambridge.

### Multi-electrode Array

Drug effect of spontaneous firings at 100 and 240 culture days *in vitro*



Epilepsy phenomenon was evoked by administration of pentylenetetrazole (PTZ) and inhibited by anti-epilepsy drug phenytoin and sodium valproate (VPA)



Experiments performed by I. Suzuki and A. Odawara at Tohoku Institute of Technology using Alpha MED Scientific Inc. MEA platform

## Conclusion

- Axol hNSCs and their progeny express neural markers at both the gene and protein level and can be cultured in 2D & 3D systems.
- Live imaging of these cells shows that they form neural networks while whole cell patch clamp and MEA recordings show the cells are electrically active.
- The properties of Axol hNSCs and hCCNs make them ideal for numerous applications including disease modeling, drug screening, toxicology studies and more.

## References

1. iPS cell technologies: significance and applications to CNS regeneration and disease. Okano & Yamanaka. Mol Brain, 2014 | 2. Induction of Pluripotent Stem Cells from Mouse Embryonic and Adult Fibroblast Cultures by Defined Factors. Takahashi & Yamanaka. Cell, 2006. | 3. Induction of Pluripotent Stem Cells from Adult Human Fibroblasts by Defined Factors. Takahashi et al., Cell, 2007 | 4. Human cerebral cortex development from pluripotent stem cells to functional excitatory synapses. Shi et al., Nature Neuroscience, 2012. | 5. A human stem cell model of early Alzheimer's disease pathology in Down Syndrome. Shi et al., Science Translational Medicine, 2012. | 6. Directed differentiation of human pluripotent stem cells to cerebral cortex neurons and neural networks. Shi et al., Nature Protocols, 2012.