

Human Dental Pulp Stem Cells

Catalog No.	Product Name	Product Quantity	Storage on Arrival	Thawing Instructions	Storage Once Thawed
ax3901	Dental Pulp Stem Cells	1 million cells/vial	Liquid nitrogen	Follow protocol	N/A
ax3902	Dental Pulp Stem Cell Culture Medium	500 mL	Aliquot and store at -20°C for up to 6 months	Thaw at 4°C or room temperature	Store at 4°C for up to 1 month
ax0044	Unlock	25 mL	Aliquot & store at -80°C for up to 6 months	Thaw at 4°C	Store at 4°C for up to 1 week
ax9009	MSC Chondrogenesis Medium	100 mL	Aliquot and store at -20°C for up to 6 months	Thaw at 4°C or room temperature	Store at 4°C for up to 1 month
ax9010	MSC Osteogenesis Medium	100 mL	Aliquot and store at -20°C for up to 6 months	Thaw at 4°C or room temperature	Store at 4°C for up to 1 month
ax9019	MSC Adipogenesis Medium for Bone Marrow Derived and Umbilical Cord Derived MSCs	100 mL	Aliquot and store at -20°C for up to 6 months	Thaw at 4°C or room temperature	Store at 4°C for up to 1 month
ax0049	Fibronectin Coating Solution	1 mL	Aliquot and store at -80°C for up to 3 months	Thaw at 4°C	Once diluted, use immediately

Lot-specific information such as donor details and passage number are stated in the Certificate of Analysis for each product.

General Recommendations:

Always count the number of viable cells after thawing

- Recommended culture vessel coating: Not required
- Recommended cell culture medium: **Dental Pulp Stem Cell Culture Medium**
- Recommended seeding density: **5,000 viable cells/cm²**
- Recommended centrifugation speed: **200 x g for 5 minutes**

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051 | US Toll Free Tel: 1-800-678-2965

Email: support@axolbio.com

Web: www.axolbio.com

Thawing & Plating:

- Transfer the vial of cells from liquid nitrogen storage with the vial buried in dry ice. Remove the vial from dry ice and transfer it immediately to a **37°C** water bath.
- Thaw the cells quickly in a **37°C** water bath until just prior to complete thawing.
- Wipe the outside of the vial with 70% ethanol.
- Gently resuspend the cells and transfer to a 15 mL sterile conical tube.
- Slowly add **10 mL** of pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium**.
- Rinse the cryovial with 1 mL of **Dental Pulp Stem Cell Culture Medium** to ensure all of the cells are transferred.
- Centrifuge the cells at **200 x g** for **5 minutes**.
- Carefully remove the supernatant and resuspend the cell pellet in **1-2 mL** of pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium**.
- Perform a cell count to determine the number of viable cells.
- Dilute the cells into the required volume of pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium**.
- Seed cells into the culture vessel at the recommended seeding density of **5,000 viable cells/cm²**.
- Incubate the cells at **37°C**, **5% CO₂** in a humidified incubator.
- After **24 h**, replace the culture medium with fresh pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium**.
- Replace the culture medium **every 1-2 days** depending on confluency.

Passaging:

- Passage when the culture reaches: **80-100% confluent**
- Recommended passaging reagent: **Unlock**
- Remove all spent medium from cell culture vessels.
- Gently rinse the surface of the cell layer once with PBS, 2 mL of PBS per 10 cm² culture surface area. Discard the PBS.
- Add 1 mL per 10 cm² of culture surface area of **cold/room temperature Unlock** passaging reagent. Evenly distribute it over the entire cell layer.
- Incubate the cells for **5 minutes** at **37°C**. Observe the cells at regular intervals for detachment from the culture vessel.
- Once the cells have detached, dilute out the passaging reagent with four volumes pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium**. For example, if 1 mL of **Unlock** is used, then add 4 mL of the medium to stop the reaction.
- Transfer the cell suspension to a sterile conical tube.
- Centrifuge the cells at **200 x g** for **5 minutes**.

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051 | US Toll Free Tel: 1-800-678-2965

Email: support@axolbio.com

Web: www.axolbio.com

- Carefully remove the supernatant and resuspend the cell pellet in **1-2 mL** of pre-warmed, **37°C, Dental Pulp Stem Cell Culture Medium**.
- Perform a cell count to determine the number of viable cells.
- Dilute the cells into the required volume of pre-warmed, **37°C, Dental Pulp Stem Cell Culture Medium**.
- Seed cells into the culture vessel at the recommended seeding density of **5,000 viable cells/cm²** or as appropriate for subsequent lineage differentiation (see below).
- Incubate the cells at **37°C, 5% CO₂** in a humidified incubator.

It is recommended that the Dental Pulp Stem Cells are not passaged more than 3 times since early passage cells perform better for subsequent lineage differentiation.

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051 | US Toll Free Tel: 1-800-678-2965

Email: support@axolbio.com

Web: www.axolbio.com

Differentiation to Adipocytes:

Recommendations:

- Recommended culture vessel coating: **Fibronectin**
- Recommended cell culture medium: **MSC Adipogenesis Medium (ax9019)**
- Recommended seeding density: **40,000-60,000 viable cells/cm²**

Coating:

- Dilute the stock **Fibronectin Coating Solution 1:100** in sterile water to make 1x working solution e.g. **100 µL** in **10 mL**.
- Coat the required number of 6-well plates with the **Fibronectin** 1x working solution. We recommend coating at a volume of **200 µL per cm²** however, please optimize for your experiments.
- Incubate the culture vessel **overnight** at **37°C** in a humidified incubator.

Coating the culture vessels with fibronectin will prevent the cells from detaching during differentiation

Adipocyte Differentiation (28 days):

- Seed cells into the **fibronectin-coated** culture vessel at the recommended seeding density of **40,000-60,000 viable cells/cm²** in pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium**.
- The cells should be maintained at **37°C**, **5% CO₂** in a humidified incubator.
- At **24 h** post-seeding, gently remove the **Dental Pulp Stem Cell Culture Medium** and rinse the wells with PBS once.
- Remove the PBS and add pre-warmed, **37°C**, **MSC Adipogenesis Medium** (2 mL per well for a 6-well plate).
- Replace the **MSC Adipogenesis Medium every 3-4 days**. Gently remove the spent cell culture medium and dispense pre-warmed, **37°C**, fresh culture medium down the side wall of the well.
- At **Day 28** after seeding, differentiation will be complete.
- Fix the cells with 4% paraformaldehyde and stain for lipid accumulation such as by using Oil Red O stain.

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051 | US Toll Free Tel: 1-800-678-2965

Email: support@axolbio.com

Web: www.axolbio.com

Differentiation to Chondrocytes:

Recommendations:

- Recommended culture vessel coating: Not required
- Recommended cell culture medium: **MSC Chondrogenesis Medium (ax9009)**
- Recommended seeding density: **1.6×10^7 cells/mL**

Chondrocyte Differentiation (21 days):

- Resuspend the cells in pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium** at a concentration of **1.6×10^7 cells/mL**.
- Seed **5 μ L** drops of cells into multi-well plates (no more than 1 drop per 0.33 cm²).
- Incubate the plates for **2 h** at **37°C**, **5% CO₂** in a humidified incubator. Ensure that the incubator is adequately humidified.
- Gently add pre-warmed, **37°C**, **MSC Chondrogenesis Medium** to the wells (2 mL per well for a 6-well plate).
- The cells should be maintained at **37°C**, **5% CO₂** in a humidified incubator.
- Gently replace the culture medium every 2-3 days with pre-warmed, **37°C**, fresh **MSC Chondrogenesis Medium**.
- At **Day 21** after seeding, differentiation will be complete.
- Fix the cells with 4% paraformaldehyde and stain for sulfated proteoglycans such as by using Alcian Blue stain.

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051 | US Toll Free Tel: 1-800-678-2965

Email: support@axolbio.com

Web: www.axolbio.com

Differentiation to Osteocytes:

Recommendations:

- Recommended culture vessel coating: **Fibronectin**
- Recommended cell culture medium: **MSC Osteogenesis Medium (ax9010)**
- Recommended seeding density: **15,000-20,000 viable cells/cm²**

Coating:

- Dilute the stock **Fibronectin Coating Solution 1:100** in sterile water to make 1x working solution e.g. **100 µL** in **10 mL**.
- Coat the required number of 6-well plates with the **Fibronectin** 1x working solution. We recommend coating at a volume of **200 µL per cm²** however, please optimize for your experiments.
- Incubate the culture vessel **overnight** at **37°C** in a humidified incubator.

Coating the culture vessels with fibronectin will prevent the cells from detaching during differentiation

Osteocyte Differentiation (21 days):

- Seed cells into the **fibronectin-coated** 6-well plate at the recommended seeding density of **15,000-20,000 viable cells/cm²** in pre-warmed, **37°C**, **Dental Pulp Stem Cell Culture Medium** (2 mL per well for a 6-well plate).
- The cells should be maintained at **37°C**, **5% CO₂** in a humidified incubator.
- At **24 h** post-seeding, gently remove the **Dental Pulp Stem Cell Culture Medium** and rinse the wells with PBS once.
- Remove the PBS and add pre-warmed, **37°C**, **MSC Osteogenesis Medium** (2 mL per well for a 6-well plate).
- Gently replace the **MSC Osteogenesis Medium every 3-4 days**, following the notes below.
- Gently remove **75%** of the spent cell culture medium and gently dispense an equivalent volume of fresh, pre-warmed, **37°C**, culture medium down the side of the well. Do not tilt the plate when removing the culture medium and do not remove all of the culture medium.

The monolayer is fragile and must be handled with extreme care or it may detach from the plate

- At **Day 21** after seeding, differentiation will be complete. Fix the cells with 4% paraformaldehyde and stain for calcium deposition such as by using Alizarin Red stain.

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051 | US Toll Free Tel: 1-800-678-2965

Email: support@axolbio.com

Web: www.axolbio.com

Usage Statement:

Our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, *in vitro* diagnostic uses, *ex vivo* or *in vivo* therapeutic uses or any type of consumption or application to humans.

Got any questions? Need help with the protocol?
Contact Axol Technical Support at support@axolbio.com
International phone **+44-1223-751-051**
US phone **+1-800-678-AXOL (2965)**

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park,
Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051 | US Toll Free Tel: 1-800-678-2965

Email: support@axolbio.com

Web: www.axolbio.com

Human Dental Pulp Stem Cell Protocol