

Human Intra-Hepatic Biliary Epithelial Cells

| Catalog No. | Product Name | Product quantity | Short-term Storage | Long-term Storage | Thawing Instructions |
|-------------|--|--------------------|--------------------|--------------------|----------------------|
| ax3765 | Intra-Hepatic Biliary Epithelial Cells | 500,000 cells/vial | Liquid nitrogen | Liquid nitrogen | See below |
| ax3766 | Intra-Hepatic Biliary Epithelial Cell Culture Medium | 500 mL | 4°C for 1 month | -20°C for 6 months | Thaw at 4°C or RT |

Lot-specific information such as donor details and passage number are stated in the Certificate of Analysis for each product.

Recommendations:

Always count the number of viable cells after thawing

- Recommended culture vessel coating: Type I Collagen
- Recommended cell culture medium: Axol Intra-Hepatic Biliary Epithelial Cell Culture Medium
- Recommended seeding density: 5,000 viable cells/cm²
- Recommended centrifugation speed: 250 x g for 5 min

Thawing & Plating:

Coat the culture vessels with Type I Collagen or use pre-coated culture vessels

- Thaw the cells quickly in a 37°C water bath until just prior to complete thawing
- Wipe the outside of the vial with 70% ethanol
- Gently resuspend the cells and transfer to a 15 mL sterile conical tube
- Slowly add 10 mL of pre-warmed **Intra-Hepatic Biliary Epithelial Cell Culture Medium**
- Rinse the cryovial with 1 mL of **Intra-Hepatic Biliary Epithelial Cell Culture Medium** to ensure all of the cells are transferred
- Centrifuge the cells at 250 x g at room temperature for 5 min
- Carefully remove the supernatant and resuspend in 1-2 mL of pre-warmed **Intra-Hepatic Biliary Epithelial Cell Culture Medium** and perform a cell count
- Dilute the cells into the required volume of pre-warmed **Intra-Hepatic Biliary Epithelial Cell Culture Medium**

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051

Email: support@axolbio.com

Web: www.axolbio.com

- Seed cells into the culture vessel (coated with type I collagen) at the recommended seeding density
- After 24 h, replace the culture medium with fresh, pre-warmed **Intra-Hepatic Biliary Epithelial Cell Culture Medium**
- Frequency of media changes: Every 2 days

Passaging:

- Passage when the culture reaches: 80% confluent
- Recommended passaging reagent: Trypsin-EDTA
- Neutralize the trypsin with **Intra-Hepatic Biliary Epithelial Cell Culture Medium** and centrifuge the cells at 250 x *g* for 5 min
- Remove the supernatant and resuspend in 1-2 mL of pre-warmed **Intra-Hepatic Biliary Epithelial Cell Culture Medium**
- Perform a cell count to determine the number of viable cells
- Dilute the cells into the required volume of pre-warmed **Intra-Hepatic Biliary Epithelial Cell Culture Medium**
- Seed cells into the culture vessel (coated with type I collagen) at the recommended seeding density

Usage Statement:

Our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, *in vitro* diagnostic uses, *ex vivo* or *in vivo* therapeutic uses or any type of consumption or application to humans.

Address: Axol Bioscience Limited, Suite 3, The Science Village, Chesterford Research Park, Little Chesterford, Cambridgeshire, CB10 1XL

Phone: +44 (0) 1223 751 051

Email: support@axolbio.com

Web: www.axolbio.com

Human Intra-Hepatic Biliary Epithelial Cell Protocol