


# **Synchronized differentiation of Axol Human Neural Stem Cells (NSCs) into Cerebral Cortical Neurons**

**Xeno-Free System**

Instruction Manual

Version 2.5

XF Protocol - 1



# Table of Contents


Product Information	3
Schematic Overview and Density Guidelines	4
Preparing Axol Plating-XF Medium	5
Preparing Axol Neural Differentiation-XF	5
Preparing Axol Neural Maintenance-XF	5
Preparing Matrix for Adherent Cell Culture Axol Sure BondXF (ax0041XF)	6
Thawing Axol NSCs for Synchronous Differentiation	7
Technical Support	9

## Product Information

Catalog no.	Product Name:	Format	Stock Concentration.	Storage on Arrival:	Thawing Instructions:	Storage Once Thawed:
<b>ax0034-125</b>	Axol Neural Differentiation-XF Medium	1 x 125 mL	1X	Aliquot and store at -80°C for up to 6 months. Keep in dark	Overnight at 4°C	Once, thawed, store aliquot at 4-8°C for up to 1 week
<b>ax0032-500</b>	Axol Neural Maintenance-XF Medium	1 x 500 mL	NA	Aliquot and store at -80°C for up to 6 months. Keep in dark	Overnight at 4°C	Once, thawed, store aliquot at 4-8°C for up to 1 week
<b>ax0041XF</b>	Axol SureBondXF	1x 1 mL	200X	4°C	N/A	Store at 4-8°C for up to 1 month
<b>ax0033</b>	Axol Plating-XF Medium	1 x 30 mL	1X	-20°C	Overnight at 4°C	Must be used immediately once thawed

# Fully-defined system to synchronously

## Protocol 1: Pure neuron population



**Density Guidelines: Protocol 1**

Catalog Number	Background	Minimum Density	Maximum Density
ax0015	Healthy	50,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>
ax0016	Healthy	50,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>
ax0111	Alzheimer's Disease (ApoE4 homozygote)	60,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>
ax0112	Alzheimer's Disease (PSEN-1 L286V)	60,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>
ax0113	Alzheimer's Disease (PSEN-1 M146L)	60,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>
ax0114	Alzheimer's Disease (PSEN-1 A246E)	60,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>
ax0115	Alzheimer's Disease (PSEN-2 N141L)	60,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>
ax0211	Huntington's Disease (CAG:45)	60,000 cells per cm <sup>2</sup>	200,000 cells per cm <sup>2</sup>

## Preparation of Plating-XF Medium

---

1. Upon receipt, store **Axol Plating-XF Medium** at or below **-20°C** protected from light. Stored at **-20°C**, media is stable for 6 months from date of manufacture.
2. When ready to use, thaw plating media overnight at **4°C** in the dark.
3. Once thawed, **Axol Plating-XF Medium** should be used immediately and **should not** be used for subsequent experiments.

## Preparation of Neural Differentiation-XF Medium

---

1. Upon receipt, aliquot and store your **Axol Neural Differentiation-XF Medium** at or below **-20°C** protected from light. Stored at **-20°C**, media is stable for 6 months from date of manufacture.
2. When ready to use, thaw an aliquot of media overnight at **4°C** in the dark.
3. A thawed, supplemented aliquot of **Axol Neural Differentiation-XF** can be stored at **4°C** for 1 week. Protect from light.

## Preparation of Neural Maintenance-XF

---

1. Upon receipt, the user should aliquot and store **Axol Neural Maintenance-XF** at or below **-20°C** protected from light. Stored at **-20°C**, media is stable for 6 months from date of manufacture.
2. When ready to use, thaw an aliquot of media overnight at **4°C** in the dark.
3. A thawed, supplemented aliquot of **Axol Neural Maintenance-XF Medium** can be stored at **4°C** for 1 week. Protect from light.

# Preparing Matrix for Adherent Cell Culture (ax0041XF)

---

1. Check the total number of viable cells on the cryovial or on the Certificate of Analysis shipped with the cells.
2. Calculate the total surface area that requires coating. This is the total number of viable cells (e.g. 2 million) / your desired plating density (see page 4 for guidelines).
3. Dilute the **Axol SureBondXF** stock solution (**200X**) in D-PBS (without calcium or magnesium) to make **1X working solution e.g. 30  $\mu$ L in 6 mL**.
4. Coat the surface of your culture vessel with the **Axol SureBondXF 1X working solution**. We recommend coating at **200  $\mu$ L 1X solution per  $\text{cm}^2$**
5. Incubate for **4 hours** at **37°C**.

**Warning:** Do not wash the vessel after coating with **Axol SureBondXF**.  
Do not allow **Axol SureBondXF** -coated culture vessels to dry.

# Thawing Axol NSCs for Synchronous Differentiation

1. Remove the cells from dry ice or liquid nitrogen storage. Immediately transfer the cells to a **37°C water bath**.
2. Quickly thaw the vial of cells by swirling it in the **37°C water bath**. Do not completely submerge the vial. Remove the vial before the last bit of ice has melted.
3. When thawed, immediately transfer the cells into a 15 mL sterile conical tube, and carefully add **10 mL** of **Axol Plating-XF Medium**.
4. Centrifuge the cells at **200 g for 5 mins**, and discard the supernatant.

**Please count cells to ensure optimal seeding density.**

5. Resuspend the cell pellet in the required amount of **Axol Plating-XF Medium**. We recommend the use of **200 µl Axol Plating-XF Medium per cm<sup>2</sup>**.
6. Quickly remove the diluted **Axol SureBondXF** coating solution from the pre-coated culture vessel before plating resuspended cells.
7. Plate the resuspended cells with **Axol Plating-XF Medium** according to density guidelines (cell line dependent) on your **Axol SureBondXF** coated culture vessel.

**Top Tip:** It is critical to promote consistent cell density, monolayer and health throughout the culture to avoid edge effects and variations in cellular maturity. After seeding, avoid disturbing the culture vessel for a minimum of 30 minutes to allow the cells to adjust to their environment.

8. Incubate the plated cells at **37°C, 5% CO<sub>2</sub>**.
9. **24 hours** after plating, replace the spent medium with fresh pre-warmed **Axol Neural Differentiation-XF Medium**.
10. After **72 hours**, replace half the volume of spent medium with fresh, **Axol Neural Maintenance-XF Medium**.

11. After **24 hours**, replace half the volume of spent medium with fresh, pre-warmed **Axol Neural Maintenance-XF Medium**. Repeat this process of re-feeding cultures with half-volumes of fresh, pre-warmed **Axol Neural Maintenance-XF Medium** every four days.

**Top Tip:** Cultures can be maintained under these conditions for over 50 days in culture!



# Technical Support

---

- Online Resources

Please visit our website at [www.axolbio.com](http://www.axolbio.com) for additional product information and *Technical Resources*, including instruction manuals, application protocols, video guides, wall charts and webinars.

- Contact Us

For more information or technical assistance, call +44 (0) 1223 497 119, or email [support@axolbio.com](mailto:support@axolbio.com). US Toll Free Tel: 1-800-678-2965 (1-800-678-AXOL), US Toll Free Fax: 1-800-861-2965 (1-800-861-AXOL).


Where we are:

[Axol Bioscience Ltd, Babraham Research Campus, Cambridge, CB22 3AT, United Kingdom.](#)

- Certificate of Analysis

The Certificate of Analysis provides detailed quality control information for each product. Certificates of Analysis are available on our website.

Go to [www.axolbio.com/certificate-of-analysis-lookup](http://www.axolbio.com/certificate-of-analysis-lookup) and search for the Certificate of Analysis with product lot number, which is printed on the cryovial label.



Here For You. We are Axol Bioscience.  
[www.axolbio.com](http://www.axolbio.com)