

KI67 Monoclonal Antibody

KI67 Monoclonal Antibody™ is a highly purified monoclonal antibody against the target protein KI67.

Product	Specification
Catalogue number	ax14001
Application	IHC-p, IF
Clonality	Monoclonal
Form	Liquid
Host	Mouse
Isotype	IgG
Shipping conditions	Blue ice
Size	60µL / 120µL
Species reactivity	Human, Mouse, Rat
Specificity	KI67
Storage buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.4
Storage conditions	-20°C
Label	Unconjugated
Formulation	1mg/mL
Clone number	Clone:8K5
Purify	Protein A purification

Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Dilution for IHC application 1:100-1:300

Dilution for IF application 1:50-200

Technical Resources

- Application Notes: www.axolbio.com/page/application-notes
- Product Page: <https://www.axolbio.com/shop/product/ki67-monoclonal-antibody-87561>
- <https://www.axolbio.com/shop/product/ki67-monoclonal-antibody-87779>